


ri-gital®

 **Riester**

The familiar way



 **Easy to use, robust and precise –**

Clinical thermometer in the well-known Riester quality. Can be used for oral, rectal and under-arm measurement of body temperature. Digital measurement is particularly accurate thanks to the highly sensitive temperature sensor.

- Maximum thermometer
- With 3-digit display, memory and signal tone
- Unbreakable and water-tight
- Measurement range: 32 – 43.9°C
- Accuracy: +/- 0.1°
- Automatic switch-off after 10 min.

- Incl. long-life button cell batteries
- In plastic storage box

ri-gital® Digital Thermometer.

No. 1850

Packet with 40 disposable probe sheaths.

No. 12703