

precisa® N

**Riester**

The familiar way



A typical **Riester quality product**. Available in metal with an aluminium housing and simple tube - plug-in connection or alternatively with a sturdy plastic housing with chromium-plated tube - screw connection. precisa® N is ideal as a mobile, secondary or reserve device for any office practice and every hospital. Also optionally available in a latex-free version (when ordering enter LF for latex-free).

- Latex bulb for pumping up the cuff
- Metal valve base
- Ergonomic spoon-shaped grip of stainless steel
- Single-tube connection
- Precision air release valve with fine adjustment and wear-free
- Microfilter protects air release valve and measuring system
- Specially tempered copper-beryllium diaphragm, with precision movement and almost non-aging
- Diaphragm pressure-loading capacity up to 600 mm Hg
- Linear, tilted aluminum scale, from 0 to 300 mm Hg (easily legible)
- Maximum error tolerance +/- 3 mm Hg
- Comprehensive selection of different sizes and versions of cuffs

	Ø 64 mm aluminium		Ø 63 mm plastic/metal	
	Single-tube	Twin-tube	Single-tube	Twin-tube
Cotton bandage cuff, adult	<b>No.1362-101</b>	<b>No.1363-103</b>	<b>No.1360-101</b>	<b>No.1361-103</b>
Cotton hook cuff, adult	<b>No.1362-104</b>	<b>No.1363-106</b>	<b>No.1360-104</b>	<b>No.1361-106</b>
Calibrated nylon velcro cuff, adult	<b>No.1362</b>	<b>No.1363</b>	<b>No.1360-107</b>	<b>No.1361-109</b>